NobelProcera™ Protocole clinique de scellement

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Conventional composite-resin luting agents and conventional silane coupling agents do not provide long-term bond strengths to high-strength ceramics.

A primer or resin cement that contains special adhesive monomers that have the ability to chemically bond to metal oxides are needed.

Blatz MB, Oppes S, Chiche GJ, Holst S, Sadan A. Influence of Cementation Technique on Fracture Strength and Leakage of Alumina All-Ceramic Crowns after Cyclic Loading.

Quintessence Int 2008;39:23-32.

The use of the MDP-containing composite resin on air abraded zirconia ceramic can be recommended as promising bonding method.

Blatz MB, Sadan A, Martin J, Lang B. In-vitro evaluation of shear bond strengths of resin to densely-sintered high-purity zirconium-oxide ceramic after long-term storage and thermocycling. J Prosthet Dent 2004;91:365-362.

Wolfart M, et al. Durability of the resin bond strength to zirconia ceramic after using different surface conditioning methods. Dent Mater. 2007;23:45-50.

Kern M, Wegner SM. Bonding to zirconia ceramic: adhesion methods and their durability. Dent Mater. 1998;14:64-71.

Atsu SS, Kilicarslan MA, Kucukesmen HC, Aka PS. Effect of zirconium-oxide ceramic surface treatments on the bond strength to adhesive resin. J Prosthet Dent. 2006;95:430-436.

Valandro LF, Ozcan M, Bottino MC, Bottino MA, Scotti R, Bona AD. Bond strength of a resin cement to high-alumina and zirconia-reinforced ceramics: the effect of surface conditioning. J Adhes Dent. 2006;8:175-181.

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Silanising sand blasted Al_2O_3 increased bond strength of conventional resin cements significantly.

Sadan A, et al. Influence of silanization on early bond strength to sandblasted densely sintered alumina. Quintessence Int. 2003;34:172-176.

Blatz MB, Sadan A, Arch G, Lang B. In vitro evaluation of long-term bonding of Procera® AllCeram alumina restorations with a modified resin luting agent. J Prosthet Dent 2003;89(4):381-387.

Ceramic cleaning methods after try-in procedures have a significant influence on the resin-ceramic bond strength. Air abrasion of contaminated zirconia ceramic is the most effective.

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